CHLORIDE

Cat. No.	Pack Name	Packaging (Content)
XSYS0008	CL 120	R1: 10 x 12 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Chloride in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Chloride is estimated in conjunction with the other electrolytes and used to support interpretation of the other electrolytes. The anion gap ($\text{Na} + \text{K} - (\text{Cl} + \text{HCO}_3)$) between these electrolytes is frequently affected in disease states. Increased levels can occur in nephritis, prostatic obstruction, eclampsia and dehydration. Decreased levels can occur with impaired gastrointestinal or renal function.

PRINCIPLE

When chloride is mixed with a solution of undissociated mercuric thiocyanate, the chloride preferentially combines with mercury forming mercuric chloride. The thiocyanate that is released then combines with ferric ions present in the solution forming strongly coloured ferric thiocyanate with an absorption maxima at 480 nm.

REAGENT COMPOSITION

R1	
Mercuric Thiocyanate	1.01 mmol/l
Ferric Nitrate	37.13 mmol/l
Nitric Acid	62.43 mmol/l
Mercuric Nitrate	105 µmol/l

REAGENT PREPARATION

Reagent R1 is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION & HANDLING

Use serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability: 7 days at 20–25 °C
7 days at 4–8 °C
at least one year at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to coulometry.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

Results are calculated automatically by the instrument.

UNIT CONVERSION

mmol/l = mEq/l

EXPECTED VALUES ¹

Serum:

Adult	98–107 mmol/l
Cord	96–104 mmol/l
Premature	95–110 mmol/l
0–30 d	98–113 mmol/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 3.7 mmol/l

Linearity: 160 mmol/l

Measuring range: 3.7–160 mmol/l

Intra-assay precision Within run (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	116.1	1.25	1.08
Sample 2	108.3	1.64	1.52

Inter-assay precision Run to run (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	117	1.47	2.75
Sample 2	94.6	1.43	1.52

COMPARISON

A comparison between XL-Systems Chloride (y) and a commercially available test (x) using 40 samples gave following results:

$y = 1.000 x + 1.0 \text{ mmol/l}$

$r = 0.927$

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 500 mg/dl.

Note:

Grossly haemolysed, lipaemic and jaundiced samples give falsely elevated results and should not be used. High levels of immunoglobulins, as found in case of multiple myeloma will interfere due to the development of turbidity.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains < 6 % methanol.

Danger

Hazard statement:

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

H371 May cause damage to organs.

Precautionary statement:

P260 Do not breathe vapours.

P280 Wear protective gloves/protective clothing/eye protection.

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 IF exposed or concerned: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units / SI units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	CLO	CLO	CLO	CLO	CLO	CLO
Test Code	49	49	49	49	49	49
Report Name	Chloride	Chloride	Chloride	Chloride	Chloride	Chloride
Unit	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	505	505	505	505	505	505
Wavelength-Secondary	660	660	660	660	660	660
Assay type	1-Point	1-Point	1-Point	1-Point	1-Point	1-Point
Curve type	Linear	Linear	Linear	Linear	Linear	Linear
M1 Start	0	0	0	0	0	0
M1 End	0	0	0	0	0	0
M2 Start	5	5	6	9	7	5
M2 End	8	8	12	15	8	8
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	1.6	1.6	1.6	1.6	1.6	1.6
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	3.7	3.7	3.7	3.7	3.7	3.7
Technical Maximum	160	160	160	160	160	160
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	0	0	0	0	0	0
Reagent Abs Max	0.3	0.3	0.3	0.3	0.3	0.3
Auto Rerun	No	No	No	No	No	No
Total Reagents	1	1	1	1	1	1
Reagent R1	CLO R1	CLO R1	CLO R1	CLO R1	CLO R1	CLO R1
Reagent R2	NA	NA	NA	NA	NA	NA
Reagent R3	NA	NA	NA	NA	NA	NA
Test Volumes						
Test	CLO	CLO	CLO	CLO	CLO	CLO
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes						
Normal	2	2	2	2	2	2
Dilution Ratio	1	1	1	1	1	1
Increase	4	4	6	4	4	4
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	2	2
Dilution Ratio	5	5	5	5	5	5
Standard volume	2	2	2	2	2	2
Reagent Volumes and Stirrer speed						
RGT-1 Volume	200	200	200	200	200	200
R1 Stirrer Speed	High	High	NA	High	High	High
RGT-2 Volume	0	0	0	0	0	0
R2 Stirrer Speed	NA	NA	NA	NA	NA	NA
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
Reference Ranges						
Test	CLO	CLO	CLO	CLO	CLO	CLO
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	98	98	98	98	98	98
Normal-Upper Limit	107	107	107	107	107	107
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	98	98	98	98	98	98
Normal-Upper Limit	107	107	107	107	107	107
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Revision Number						
Revision	<A-100- CLO-1 20.08.2013>	<A-200- CLO-1 20.08.2013>	<A-300/600- CLO-1 20.08.2013>	<A-640- CLO-1 20.08.2013>	<A-1000- CLO-1 20.08.2013>	<A-180- CLO-1 12.12.2013>

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1. Tietz Textbook of Clinical Chemistry and Molecular diagnostics. Burtis, C.A., Ashwood, E. R., Bruns, D.E.; 5th edition, WB Saunders Company, 2012.
2. Zall, D.M. Fischer, D. Garner, D.O. Anal. Chem. 28, 1665, 1956.
3. Harper, R.J. (Ed.) Clinical Chemistry: Principles and Techniques (2nd Ed.) Harper and Row 1974, P. 718 - 719.

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